

Where to place markers


Correct positioning of the markers on your subject is key in achieving good motion capture. Although you could place the markers anywhere, bear in mind that you are usually aiming to capture the movement of the underlying skeleton.

As such, markers should be placed at specific and easily identified anatomical locations.

Where to attach markers


The figure below shows an example of where to position markers for a full body capture

Front and rear views for the example marker placement. Note the asymmetry in upper arm, forearm, thigh and shin marker.



The names of the markers are called labels. These labels are listed in a file called the Marker Set and are associated with the markers in the Workstation application. Marker sets are described in greater detail in the Reference section.

Hand:



RIGHT HAND